

NEW




ACP-1500

Auto Chart Projector

ACP-1500

Projection distance	1.5~7 meter
Projection magnification	30×(at 5m)
Projection size	330×270mm (at 5m)
Chart Switching speed	0.01~0.02S
Mask number	1 open, 5 horizontal lines, 8 vertical lines, 21 single letters,1 red/green
Mask Switching speed	One mask per 0.01~0.02S
Program	2 sets programs, each program contains up to 30 steps
Lamp	12V、50W (halogen lamp) 12V、10W(LED)
Power saving function	Will automatically standby in 10 min
Power source	AC 220V, 50Hz or 110V,60Hz
Power consumption	80 W
Dimensions	320mm ×218mm×240mm
Net weight	7.2 KG
Standard accessories	Remote control, polarized metal screen, halogen lamp, polarized glasses,fuses(2), battery(2)
Optional accessories	Floor stand,Wall stand



Polarized glasses

ACP – 1500B	E 400 120	W 300 90	END 200 60 SLH 160 45	HBDVA 100 30 PHTCK 80 24	VLNEB 70 21 DAOFP 60 18 EGDNS 60 15	FZBOT 40 12 OFLCV 30 9 APECH 25 7,5	
TZVFC 20 6 EHPNL 15 4,5 VEOTB 10 3	EVOTL 20 6 TEGAK 20 6 HSFZN 20 6	WEM 200 60 M3W 100 30	WEM3W 80 24 EM3W3 70 21 M3WEM 60 18	EW3MW 50 15 W3ME3 40 12 MEW3M 30 9	3WEMW 25 7,5 WEM3E 20 6 EM3W3 15 4,5	852 200 60 643 160 45	29387 100 30 75462 80 24 38549 70 21
69832 60 18 42576 50 15 93754 40 12	38297 30 9 53948 25 7,5 82573 20 6	200 60	100 30 80 24	60 18 40 12 30 9			
i-Optik ACP – 1500C	K 0.05	DEC 0.1 HON 0.16	TKPEB 0.2 ARFSH 0.3 CNDTZ 0.4	KOZFT 0.5 PBHAD 0.6 RZCNB 0.7	VSHLE 0.8 APFCB 0.9 RDZTG 1.0	HKNST 1.2 OVEFG 1.5 DAGZR 2.0	
28563 0.6 93427 0.7 56984 0.8	94527 1.0 68349 1.2 59283 1.5	0.1	0.2 0.3	0.4 0.5 0.6	0.7 0.8 1.0		