$\epsilon$ 











## **General Features**

5-inch touch screen 11 testing parameters 100 samples/hour throughput Dual wavelength reflectance test method Auto & manual test mode Thermal printer, auto & manual print Support LIS and barcode Less and easy maintenance

## **Technical Specifications**

Principle: Dual wavelength reflectance method
Wavelength: 520nm, 640nm, 660nm
Parameters: Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein,
Urobilinogen, Nitrite, Leukocyte, Vitamin C

Throughput: 100 samples/hour Sample type: Urine Test mode: Auto, manual Waiting time: 15-40s Display: 5-inch touch screen Memory: 2,000 results storage

Memory: 2,000 results storage
Language: English
Strip: Original strip, 100 strips per bottle
Report format: Semi-quantitative and qualitative
Coincidence rate: ≥90%
Power supply: AC 100-240V 50/60Hz
Working temperature: 10-30°C

Atmospheric pressure: 86.0kPa-106.0kPa

Relative humidity: ≤70% Size: 360mm×225mm×160mm

Net weight: 2.5kg



