










1, 3, 6 Channels ECG Machine Comparison

	One Channel		Three Channels		Six Channels
	ECG-101	ECG-101G	ECG-300	ECG-300G	ECG-6010
Image					
Screen	Smaller LCD Single Waveform Display	Larger LCD Three Waveforms Display	Smaller LCD Three Waveforms Display	5" LCD Twelve Waveforms Display	5.7" TFT LCD Twelve Waveforms Display
Operation	By Key				By Key and Shuttle
Storage	No	Optional 2G TF Card	20	Optional 2G TF Card	Optional 2G SD Card
Language	Russian, English, Spanish, French, German, Italian				
Interface	RS232				USB, LAN
Paper	50mm X 30m		63mm X 30m		112mm X 30m (Roll) 112mm X 140mm X 20m (Fold)
Interpret Function	No	Yes			Advanced
Dimension	310mm X 230mm X 70mm				290mm X 245mm X 99mm
G.W.	2.7 Kg		2.8 Kg		4.6 Kg

12 Channels ECG Machine Comparison

	ECG-1200	ECG-1210	ECG-1215	ECG-1230
Image				
Screen	5.7" LCD Resolution: 320 X 240	5.7" TFT LCD (Foldable) Resolution: 320 X 240	7" TFT LCD (Foldable) Resolution: 800 X 480	12.1" TFT Color LCD Resolution: 800 X 600
Screensaver	No	No	Yes	Yes
Operation	By Key	By Key and Shuttle		
Storage	>800	Optional 2G SD Card		
Language	Russian, English, Spanish, French, German, Italian			
Interface	RS232	USB,LAN		
Recording Paper	210mm×30m (Roll)	210mm×140mm×20m (Fold)		
Interpret Function	Advanced			
Dimension	400mm×325mm×84.5mm	360mm×325mm×86mm		360mm×325mm×380mm
G.W.	5Kg	6.6 Kg	6.6 Kg	8.6 Kg