

Dry Oven

DHG series

- ◆ Microprocessor control
- ◆ With double glass window
- ◆ Digital display setting and actual temperature
- ◆ Digital display running timer
- ◆ Stainless-steel interior, chamber corners are semicircle

Safety Mechanism:

- ◆ PID self-tuning, over temperature alarm
- ◆ Automatic alarm after operation

Scope of application:

- ◆ This product is suitable for drying, baking, melting wax sterilization and curing in industrial and mining enterprises, medicine and health, and scientific research units.



DHG-9053A

DHG-9423A

DHG-9240A



Model	DHG-9023A	DHG-9053A	DHG-9123A	DHG-9423A
Capacity	25L	50L	105L	420L
Working size(mm)	300x300x275	415x375x345	550x370x520	600x550x1300
Power	500W	1000W	1500W	4000W
Temperature range	50~200°C			
Temperature resolution	0.1°C			
Temperature fluctuation	±1°C			
Timer	0-9999min			
Voltage	110V 60Hz / 220V 50Hz			
Dimension(mm)	570x530x460	690x610x530	720x700x720	870x1000x1500
Package(mm)	670x610x530	770x670x660	920x730x820	850x920x1750
G.W./N.W	41/29kg	56/40kg	78/58kg	210/170kg

Model	DHG-9030A	DHG-9070A	DHG-9140A	DHG-9240A	DHG-9624A
Capacity	40L	70L	140L	225L	640L
Working size(mm)	350x350x350	450x350x450	550x450x550	500x600x750	800x800x1000
Power	500W	1000W	1500W	2000W	5000W
Temperature range	50~300°C				
Temperature resolution	0.1°C				
Temperature fluctuation	±1°C				
Timer	0-9999min				
Voltage	110V 60Hz / 220V 50Hz				
Dimension(mm)	620x630x520	620x710x620	720x780x750	770x840x980	940x1090x1380
Package(mm)	680x720x610	720x750x740	830x860x870	900x870x1130	1000x1190x1510
G.W./N.W	52/38kg	62/48kg	86/64kg	111/89kg	240/200kg