

## **DKS Series Water Bath**

## [Uses Outlined]

DKS series electric heating constant temperature water-bath water pot can supply medical health, medical colleges and check the pathological serum and industrial and mining enterprises, scientific research units do precision constant temperature and auxiliary heat It adopts stainless steel core and movable top lid. It is precise in temperature control and is well equipped with digital temperature display and chasing alert system.



## [DKS Series Water bath]

Туре	DKS-11	DKS-12	DKS-14	DKS-16	DKS-24	DKS-26	DKS-28
Input Power	220V 50Hz						
Consumption Power	300W	500W	1000W	1500W	1000W	1500W	2000W
Temperature range	RT+5°C ~ 99°C						
Temperature fluctuation	±1℃						
Follow-up alarm	+2°C						
Working Size(mm)	150×150×110	300×150×110	600×150×110	900×150×110	300×300×110	450×300×110	600×300×110
Overall Size(mm)	300×192×195	450×192×195	750×192×195	1050×192×195	450×345×195	600×345×195	750×345×195