

Vet Electrocardiograph



ECG-101G

Features

- Digital isolation technology minimized signal interference;
- Wide screen LCD with 3 lead or 1 lead waveforms display;
- 5 leads ECG cable specially for veterinary using;
- Build-in rechargeable battery, AC/DC supply;
- Power off automatically after long time non-operation;
- The width of baseline is adjustable, convenient to observe;
- 7 leads recording by high resolution and precision thermal array printer;
- Internal memory supporting ECG data and waveforms store and recall.



ECG Cable



ECG Cable 4.6m
(For horse, ect.)



ECG Electrode

Functional upgrade are available.
For detail information, please contact us.

