

# ECG-1200

Digital 12-Channel ECG



# ECG-1200

## Digital 12-Channel ECG

### Features:

- 5.7" high brightness LCD display 12 waveforms;
- R-R interval analysis, detailed and brief different analysis reports;
- Various work modes satisfy different requirements;
- Economy mode available help to reduce printing cost;
- High ability Li-ion battery up to 2 hours continues working;
- Multi-language options;
- RS232 interface communicate with PC;
- Internal memory supporting the maximum 2GB data storage.

### Specifications:

|                         |  |
|-------------------------|--|
| Lead                    | Standard 12 leads, simultaneous acquisition                          |
| Input Circuit           | Floating circuit input   |
| Input Impedance         | $\geq 50M\Omega$   |
| Frequency Response      | 0.05~150Hz   |
| Patient Leakage Current | $< 10\mu A$  |
| Input Circuit Current   | $\leq 0.1\mu A$  |
| Time Constant           | $\geq 3.2s$  |
| Sensitivity             | 2, 5, 10, 20mm/mV  |
| Paper Speed             | 5, 6.25, 10, 12.5, 25, 50mm/s  |
| Operation Mode          | Automatic, manual and rhythm   |
| Power Requirements      | AC: 100~240V, 50/60Hz, 55VA<br>DC: 14.8V rechargeable Li-ion battery |
| Safety Level            | IEC Class I, type CF   |
| Recording Paper         | 210mm $\times$ 140mm $\times$ 20m, fan-fold paper (A4 size)          |
| Dimension of Packing    | 500mm $\times$ 420mm $\times$ 260mm (L $\times$ W $\times$ H)        |
| Net Weight              | 5Kg  |



Inner Handhold



Standard Accessories



### Pioway Medical Lab Equipment Co., Ltd.

Add: Room 515A, Xinwanli Building, No. 9 Baixia Road, Nanjing, China  
Tel: +86-25-84677741 Fax: +86-25-84677742  
<http://www.pioway.com> E-mail: [export@pioway.com](mailto:export@pioway.com)