


# ECG-1215

Digital 12-channel ECG



# ECG-1215

## Digital 12-channel ECG



⦿ Detailed Analysis Report

⦿ LAN Port

⦿ Foldable Color LCD

⦿ ECG Information System



### Features:

- The foldable monitor(0~80 degree) with 7 inches color LCD;
- 12+1 waveforms display with grid line;
- Real-time waveform freezing, replaying, printing;
- Various analysis reports satisfy different requirement;
- Optional SD card, unlimited data storage;
- LAN port for real-time data transmission;
- Multi-language options;
- Long lasting Li-ion battery up to 2 hours continues working;
- Easy upgrading, keeps you updated.

### Specification:

Lead .....	Standard 12 leads,simultaneous acquisition
Input Circuit .....	Floating circuit input
Input Impedance .....	>50MΩ
Frequency Response .....	0.05~150Hz
CMRR .....	>100dB
Patient Leakage Current .....	<10μA
Input Circuit Current .....	≤ 0.1μA
Time Constant .....	≥ 3.2s
Sensitivity .....	2.5,5,10,20mm/mV

Paper Speed .....	5,10,25,50mm/s
Operation Mode .....	Automatic, manual and rhythm
Safety Level .....	IEC60601-1:2005 Class I Type CF
Power Requirements .....	AC:100~240V,50/60Hz,50VA DC:14.8V rechargeable Li-ion battery
Recording Paper .....	210mm×140mm×20m, fan-fold paper (A4 size)
Safety Level .....	IEC Class I, type CF
Dimension .....	360mm×325mm×86mm(L×W×H)
Weight .....	About 8.6kg

### Pioway Medical Lab Equipment Co., Ltd.

Add: Room 615A, Xinwanli Building, No. 9 Baixia Road, Nanjing, China  
 Tel: +86-25-84677741 Fax: +86-25-84677742  
<http://www.pioway.com> E-mail: [export@pioway.com](mailto:export@pioway.com)