




# ECG-6010

Digital 6-channel ECG

# ECG-6010

## Digital 6-channel ECG



Optional SD Card



Lan Port



Foldable Handhold



ECG Information System



### Features:

- Digital isolation technology and signal processing ensured the signal reliability;
- 5.7 inches high-brightness LCD;
- Detailed analysis report;
- Economy mode available help to reduce printing cost;
- Build-in rechargeable Li-ion battery;
- Optional SD card, unlimited data storage;
- LAN port for real-time data transmission;
- Multi-language options ;
- Easy upgrading keeps you updated;
- Delightful operating experience with durable silicon button, unique shuttle design.

### Specifications:

|                    |  |
|--------------------|--|
| Lead               | Standard 12 leads simultaneous acquisition                 |
| Input Impedance    | >50MΩ  |
| Frequency Response | 0.05~150Hz   |
| CMRR               | >100dB   |
| Sensitivity        | 2.5, 5, 10, 20mm/mV  |
| Power Supply       | AC:100~240V, 50/60Hz,42VA<br>DC:14.8V rechargeable battery |

|                 |  |
|-----------------|--|
| Paper Speed     | 5, 6.25, 10, 12.5, 25, 50mm/s                    |
| Operation Mode  | Automatic, manual and rhythm                     |
| Recording Paper | 112mm×27m(roll)<br>or 112mm×140mm×20m(fan- fold) |
| Safety Level    | IEC60601-1:2005 Class I Type CF                  |
| Dimension       | 290 mm×245 mm×99mm(L× W×H)                       |
| Weight          | 4.6 Kg   |