

# FM-200

Fetal Doppler



# FM-200

## Fetal Doppler



Headphone Available



Volume Adjustable



Operation Keys



Li-ion Battery



### Features:

- Signal strength indicating function for locate probe easily
- Physiological and system alarm function
- Volume is adjustable and optional for headphone
- Belt clamp for take away conveniently
- Three work mode satisfy the requirements of users
- Automatically power-off while signal or f probe failure for 2 minutes
- Portable and convenient for pregnancy

### Specification:

Ultrasonic Work Frequency: 2MHz  
Heart Rate Range: 50bpm~240bpm  
Resolution: 1bpm  
Alarm Limit: Adjustable lower limit range:50bpm~120 bpm  
Adjustable upper limit range:160bpm~240 bpm