



Lead of the premier laboratory  
**Laboratory preferred**  
 Lab Instruments Updating Products

**GK Series Dry heat sterilizer/GRX series Hot air sterilizer.**

**【Dry heat sterilizer】**

GK series Dry heat sterilizer/GRX series Hot air sterilizer. Multiply Band Programmed Controller With thermostat timing

**【Uses outlined】**

- ◎ Precision: Gk-9000 Series Dry Heat Sterilizers Using High-precision Temperature Control Liquid Crystal Display Device, A Very Humane Touch Screen Operation, Showing More Clearly Intuitive.
- ◎ Monitoring :the Device Uses Advanced Communications Technology To Support Rs 485 Or Rs232 Communication Interface, To Facilitate The Formation Of Small And Medium-sized Computer-controlled System.
- ◎ Prompt: The Dry Heat Sterilizers Advanced Far-infrared Heating Technology, The Rapid Warming, Forced Convection, Hot Air Dry After A Heat Objects Directly, Drying, Sterilization Time Was Significantly Shorter, Shortest Time To Reach The Maximum Temperature Of Only 10 Minutes.
- ◎ Safety: The Dry Heat Sterilizers Intelligent Temperature Control System, An Increase Of Independent Over- Temperature Protection, Temperature Control Instrument To Prevent The Situation Getting Out Of Control, To Protect The Disinfection And Sterilization Products And Equipment To Protect The Security Role Set Up An Independent High-power Sound, Light Alarm Device, Such As A Major Failure Of The Equipment, Can Quickly Inform The User In A Timely Manner. Over-temperature Protection, Over Temperature Limit That Is Automatically Interrupted, Tracking Alarm Function, Stirring Hot Air Circulation To Ensure That The Personnel, Equipment Safety.
- ◎ Time:to Increase Significantly The Number Of Independent Time Control Device. To Set The Time Limit, Time Of Arrival Of Automatic Shutdown, The Operators Do Not Need To Wait Next To The Equipment.
- ◎ Convenient:the Fast Sterilizer Control Equipment From Time To Time, Highly Adjustable Shelf, Stainless Steel Liner, Rounded Corners Are Easy To Clean, There Are Convenient To The Operator.



GK-9000 Series Dry heat sterilizer



GRX-9003 Series Hot air sterilizer

**【GK-9000 Series Dry heat sterilizer Technical parameters】**

Type	GK9040	GK9070	GK9140	GK9240
Main voltage	220V 50Hz			
Temperature fluctuation	50°C—250°C			
Follow-up alarm	±1°C			
Consumption power	+5°C			
Dimensions of working chamber	2000W	3000W	4000W	4500W
Working size	350×350×350	450×350×450	550×450×550	600×500×750

**【GRX-9003 Series Hot air sterilizer Technical parameters】**

Type	GRX-9013AS	GRX-9023AS	GRX-9053AS	GRX-9123AS	GRX-9203AS
Main voltage	220V 50Hz				
Temperature fluctuation	50°C—200°C				
Follow-up alarm	±1°C				
Consumption power	+5°C				
Dimensions of working chamber	500W	500W	800W	1500W	2000W
Working size	250×250×250	300×300×275	415×370×345	554×370×523	600×550×600