

Bacti-cinerator Sterilizer

IS800-A



Description:

- ◆ IS800-A Bacti-cinerator Sterilizer is used in degassing for inoculating loops and needles, instead of alcohol lamps. Due to its convenience, sterilizer is widely used in the environment of biological safety cabinet, the exhaust fan, purification table and mobile van.

Application:

Widely used in biosafety cabinets, clean bench, near exhaust fan, mobile vehicles, even outdoor and other other severe environment.

Features:

- ◆ Safety, convenient. Totally take the place of alcohol lamps.
- ◆ Only 5-7 seconds are needed, because of its highest temperature of 825°C.
- ◆ Smart, light, easy to handle.
- ◆ Can be used in anaerobic chamber.
- ◆ Podzolic organics prevent pollutions effectively.
- ◆ Smart technique makes it possible to use for a long time and much safer.
- ◆ 6-month warranty.

Parameters:

Project	Parameters
Max. temp. in central zone	825°C±25°C
Length of heating zone	140mm
Ambient temperature	5°C~40°C
Voltage	AC220V or AC110V, 50/60Hz
Dimension(mm)	120×155×180

Project	Parameters
Max. outer diameter of sterilization goods	Φ15mm
Heating time	15min
Fuse	250V 3A Φ5×20
Power	150W
Weight(kg)	1.0