

## ELECTRONIC BALANCE











The visual appearance of the product is blent into the simple concept, stylish and generous. Specially makes the fine pouring cast- aluminium outer covering, so as to enhance the balance of the anti-static anti-jamming capability. The panel with two /five operational button, Sliding glass windscreen provide big space, insightful vision and operate easily.

MODEL	JA103H	JA203H	JA303H	JA203	JA303	JA503	JA603
Capacity (g)	110	210	310	200	300	500	600
Price(USD)							
Sensor	Load Cell			Magnet			
Resolution (g)	0.001						
Min weighing (g)	0.004						
Stable Time	≤3 S						
Display	LED (Green)			LCD (white back light with black font)			
Pan Size	Φ 80 mm/ Φ 90mm						
Operation temp	5-25℃						
Repeat ability	±0.002g						
Liner	±0.003g						
Draft shield size	240mmx190mmx265mm						
Cal.weight	100g	100g 200g		200g			
Option	RS232C & RS485、 Printer						
Packing Size	Inside:415×285×430mm			Outside:475×315×485mm			
G.weight	7.5kg			10kg			