

The visual appearance of the product is blent into the simple concept, stylish and generous Sliding glass wind screen provide big space, insightful vision and operate easily. PL series is the dual range.

MODEL	JA103P	JA203P	JA303P	JA103PL
Capacity (g)	110	210	310	100/600 dual range
Price(USD)				
Sensor	Load Cell			Load Cell
Resolution (g)	0.001			0.001/0.01
Min weighing (g)	0.004			0.004/0.04
Stable Time	≤3 S			
Display	LCD (white back light with black font)			
Pan Size	Ф 90mm			Φ 90mm/ Φ110mm
Operation temp	5-25℃			
Repeat ability	±0.002g			±2d
Liner	±0.003g			±3d
Draft shield size	177mmx173mmx205mm			
Cal.weight	100g	200g		200g
Option	RS232C & RS485、 Printer、Dual display、USB			
Packing Size	Inside:350×265×345mm Outside:740×550×390mm			
G.weight	4.0 /16kg			

ANALYTICAL BALANCE



















