


[JYF-10B 壁挂式牙科 X 线机]

DENTAL X-RAY UNIT

CE
0434



性能特点

- 图像清晰，高效率，采用国际先进技术
- 低辐射，漏射线量仅为国家规定的百分之一
- 轻触按键，微电脑控制，只需一次按键即可迅速准确地确定曝光参数，操作简单
- 明室冲洗，一分钟成像，立即用于诊断，最大程度方便医生诊断
- 长达两米的吊架结构，即节省空间又方便医生进行各方位的拍摄
- 可用于口腔数字成像系统

Features

- Use the international advanced technology,high-efficiency,clear imaging
- Low leak radial,is only 1% of national standard
- Light touched button, controlled with microcomputer,only press one button to select the exposure parameter,fast and accurately
- Developing in light room,image in one minute,the doctor can diagnose conveniently for maximum
- The hanger body with 2m can save space and shoot conveniently in all direction
- It is for the use of the digital imaging system of mouth



技术参数

电源条件:	AC220V±10%, 50Hz, 1kVA
管电压:	60KV
管电流:	8mA
焦点尺寸:	0.8mm
总滤过:	2.5mmAL
曝光时间:	0.2~4 秒
泄漏辐射:	1 米外≤ 0.002mGy/h (国家标准为 0.25mGy/h)
包装尺寸:	管电流 0.5mA 焦点尺寸 0.8mm

Specifications

Power supply:	AC220V±10%,50HZ,1KVA
Tube voltage:	60KV
Tube current:	8mA
Focus size:	1.5mm
Total filtration:	2.5mmAL
Exposure time:	0.2~4s
Leak radiation:	Outside one meter ≤ 0.002mGy/h (National standard: 0.25mGy/h)
Optional:	Tube current 0.5mA Focus size 0.8mm