

PowerShocker LGT-2500 Series

The First & Unique Radial Shock Wave Therapy Device from China



Features

- Full colour 7" LCD touch-screen interface
- Large storage drawer for keeping clinical supplies
- Multilingual switch
- Heavy duty steel shell



Detachable connection socket



Storage tray for all accessories



7" LCD Touch Screen

		
Model Code	PowerShocker LGT-2500A	PowerShocker LGT-2500B
Handpiece	1 pc	2 pcs
Item Dimension	54×44.5×107.5 cm	54×44.5×107.5 cm
Item Weight	60.0 kg	65.0 kg
Input Voltage	110V/220V AC 50~60Hz	110V/220V AC 50~60Hz
Power Consumption	700 VA	700 VA
Pressure range	1.0~5.0 bar (step: 0.1 bar)	1.0~5.0 bar (step: 0.1 bar)
Frequency	1~15 Hz (step: 1 Hz)	1~15 Hz (step: 1 Hz)
Shocks	1~3000 (step: 100 shocks)	1~3000 (step: 100 shocks)
Compressor	built-in oil	built-in oil
Standard transmitters	1×15mm Radial 1×20mm Radial	2×15mm Radial 2×20mm Radial
Optional transmitters	6mm acupuncture, 15mm focus, 15mm deep	

The Supplied Transmitters

					
Diameter	D=15mm radial	D=20mm radial	D=6mm	D=15 mm focus	D=15mm deep
Indications	"Radial shock wave therapy for all tendinopathies"	"Myofascial trigger therapy, muscle- and connective tissue, dissolution of fibrosis"	shock wave acupuncture	fascial muscles neck, surface pain areas	"Deep-seated diseases, trigger, chronic diseases, other muscular problems"
Penetration depth	0-40mm	0-50mm	0-40mm	0-20mm	0-60mm
"Max. energy intensity mJ/mm ² "	0.40 mJ/mm ²	0.48 mJ/mm ²	0.32 mJ/mm ²	0.16 mJ/mm ²	0.63 mJ/mm ²

RSWT Handpiece for PowerShocker LGT-2500



Applications

- Myofascial trigger points
- Disorders of tendon insertions
- Activation and tightening of muscle and connective tissue
- Acupuncture shock wave therapy
- Cellulite treatment
- Elimination of lymphatic congestion
- Regulation of muscle tone
- Medical massage

