





Pressure Activated Safety Lancets

SIMPLE and **SAFETY** of use is what safety lancets are all about. The complete range comprises 4 versions, with a variety of needle, blade gauge and penetration depths to meet all your capillary blood sampling requirements.








Please read instructions carefully before use. Do not use if protective cap has been previously removed!


New Design - Great Feature Of Only Single Use



Carefully pull off the protective cap.



Push the firmly onto the chosen site to fire.



Dispose of the in a suitable sharps container.

Reorder#	Type	Color	Gauge/ Depth	Package
ZWII-022818		Violet	28g/1.8mm	100/box 2000/case
ZWII-022618		Yellow	26g/1.8mm	
ZWII-022324	= =====================================	Green	23g/2.4mm	
ZWII-022118		Pink	21g/1.8mm	
ZWII-022124	3 4 7	Orange	21g/2.4mm	
ZWII-021818	Neonatal	Green	18g/1.8mm	