

PW-302H

Pulse Oximeter

FEATURES

- Common AA alkaline or rechargeable batteries to support 7 continuous working hours
- Compact, small, light, easy for carrying and handling
- Large bright LED display of measurements, easy to read

Suitable for adult, pediatric, infant and neonate patients.

Adjustable audible and visual alarm

Applicable to animals selectively

■ Battery capacity indicator

Automatic power-off



Pulse Oximeter

Specifications

Patient Range

Adult, pediatric, infant, neonate patients

Sp₀₂

Measurement range: 35%~99%

Accuracy: ± 2%(70%~100%),0%~69% unspecified

Resolution:1%

Pulse Rate

Measurement range:30bpm ~ 250 bpm

Accuracy: ±2 bpm Resolution: 1bpm

Alarm

mode:2-level audio and visual alarm
Setup:User-adjustable high and low limits

Display Mode

SpO2 display :2-digit and 7-segment red LED nixie tubes PR display :3-digit and 7-segment red LED nixie tubes Pulse strength/ Battery capacity display:10-segment andbar-graph red LED nixie tube

Battery indicator: Red and Dual-color LED

Dimension

63 (W) × 120 (H) × 32 (D) mm

Net Weight

130g (not include sensor and batteries)

SELECTABLE SENSORS

pediatric soft tip SpO2 sensor

Neonate/infant wrap SpO2 sensor

Ear clip SpO₂ sensor



Environment Requirements

Temperature:

Working:0°C ~ 45°C

Transport and storage: -20°C ~ 60°C

Humidity:

Working: 15% ~ 95% (non-condensing)

Transport and storage: 10% ~ 95%

(non-condensing)

Altitude:

Working:86 kPa ~ 106 kPa

Transport and storage:50kPa ~ 106 kPa

Power Requirements

Battery type:common 1.5V AA alkaline or rechargeable batteries

Run time: 7 hours