



PW-9618

Digital Ultrasound Scanner



PW-9618

Digital Ultrasound Scanner



► Main features:

- Full digital beam-forming and dynamic aperture technology ensure excellent image quality
- 80 element, 20 channels, dual sockets
- Scanning angle 8-16 level adjustable
- Gamma correction, 8-segment GTC
- Automatic report generation(normal/OB)
- Abundant software packages for full clinic applications such as abdominal, Obstetric, gynecology, cardiology, urology, small parts, etc.
- Ergonomic design, user-friendly

► Specifications:

Scanning: electronic convex, linear and micro-convex
Display: B, B/B, B/M, M, 4B
Image conversion: up/down, right/left, black/white
Notation: real-time clock, date, ID, gender, age, focus, hospital, doctor
Scanning depth: 220mm(max)
Gray Scale: 256
Display: 10" SVGA
Video output: PAL-D
Body marks: 40
Cine-loop: 409 frames(max)
Permanent image storage: 32 frames
Image magnification: x1.0, x1.2, x1.5, x2.0
Picture in picture: 1:2
Preset color code: 0-7
Measurement of distance, circumference, area, volume, HR, EF, EDD
USB2.0 interface

► Standard probe

3.5MHz multi-frequency convex



3.5MHz Convex

► Other options

Video printer, trolley cart, acupuncture bracket & needle



► Optional probe

7.5MHz multi-frequency linear
7.5MHz multi-frequency endorectal linear
6.5MHz multi-frequency transvaginal convex
5.0MHz micro-convex



7.5 Endorectal



5.0MHz Micro-Convex



6.5 Transvaginal



7.5MHz Linear

※ (Specification subject to change without prior notice)



Pioway Medical Lab Equipment Co., Ltd.

Add: Room 615A, Xinwanli Building, No. 9 Baixia Road, Nanjing, China
Tel: +86-25-84677741 Fax: +86-25-84677742
www.pioway.com export@pioway.com

