



RE902-C-v

Anesthesia Machine

NEW



RE902-C-v Anesthesia machine

● Feature

CE Certificated

Module rack with polyester medical material, high strength and light weight.

Superior thin spray technology, delicate color, beautiful and durable.

High precious vaporizer, Halothane, Enflurane, Isoflurane and Sevoflurane could be choosed freely

Main unit are only 2Kg(w/o vaporizer), light weight and easy to move

With a dedicated handle, could move by hand-pull.

Two gas source:O₂ suitable for more patient type.

● Technical Specification

Closed、 Semi-closed、 Semi-open

Application:2~100Kg animal

Gas source: O₂ 0.25Mpa~0.65Mpa

Two-tube (O₂) : Scope, 0—1 0L/min. With linkage device,O₂ concentration not low than 25%.

O₂ flush: more than 25L/min

Respiratory mode: Manual、 Oxygen

TV: 50~2000mL

I:E: according to the doctor operation

BPM: according to the doctor operation

APL Valve pressure: 0.5-6KPa

Monitoring: Airway pressure、 respiratory mode

Three Gauges: two for O₂ , one for airway pressure

● Environment requirement

Temperature: Operation: 18-36°C Storage: -20-55°C

Humidifier: Operation: 45%-75% Storage: ≤93%

Atmospheric pressure: Operation: 0.95~1.05 bar

Storage: 500hpa~1060hpa,

Altitude: Operation: 500 to 800 mmHg (3565 to -440m)

Storage: 375 to 800 mmHg (5860 to -440m)

● Size: 380mm(w)×360mm(H)×180mm(L)

Weight: 8.5Kg



Area Agent: