



## Blood Bag Tube Sealer

### Application:

The blood bag tube sealer is a compact equipment to seal the blood bags, infusion bags, urine bags by high frequency sealing system. It's widely used in blood-bank, hospital, biological industry, etc.

### Features:

- \* No hemolysis.
- \* With suction feet to avoid moving, can be used in blood collecting vehicle.

### Technical Parameters:

Sealing time	0.8~1.2s
Sealing type	Automatic
Diameter of tube	Ø2-6mm (tube wall thickness≤1mm)
Power Supply	AC110V/220V±10%, 50/60HZ
External Size (W*D*H)	320*180*130mm
Package Size (W*D*H)	410*210*210mm
Net/Gross Weight	5kg/6kg

