





Excellent Performance

- Compatible with U, V or flat-bottom plates
- Up to 50 wash protocols programmable
- Positive displacement pump with low noise



Easy to Use

- 8-way and 12-way wash heads interchangeable
- Unattached absorbing head and dispensing head
- Bottom washing function



Good Cleaning Effect

- Minimal residual volume
- Bottom washing function
- Overfill washing and overflow protection



Reliable Spare Parts Quality

- Well-know spare parts supplier
- High performance PCB board
- Individual quality control management for each single component

