

PRECISION BALANCE



XY-2C

Easy to read LCD display, one touch tare up to full capacity, high resolution , quick readable weighting,easy to operate, widely used in high precise industrial applications, laboratory, university.

MODEL	XY2002C	XY3002C	XY5002C	XY6002C	XY10002C	XY20002C	XY30002C
Capacity (g)	210	310	510	610	1100	2100	3100
Price(USD)							
Resolution (g)					0.01		
Min weighing (g)					0.04		
Stable Time					≤ 2 S		
Pan Size					Φ 130 mm		
Operation temp					5-25°C		
Repeat ability					0.02g		
Liner					0.03g		
Cal.weight	200g	200g	200g	200g	500g	1000g	2000g
Option					RS232C & RS485、 Printer		
Packing Size			Inside:310×230×130mm		Outside:660×470×320mm		
G.weight					2.0/22kg		

* The calibration weight is optional parts when the capacity ≥2000g or resolution >0.01g.