

# NAVIGATION TECHNOLOGY TO BENEFIT MANKIND

## Water Distiller

### 【Uses Outlined】

The stainless steel electric automatic distillation Water distiller is specially designed for medicine, plating, lab and sanctuary use. Using tap water, it is safe, stable, capacious, properly structured, and easy to use.

### 【YNZD Water distiller】

Evaporating pot is made of quality stainless steel with all the Argon welding surfaces safe and smooth.

### 【YNSD Automatic cut-off water distiller】

- ◎ This product adopts quality stainless steel and is made out of advanced technology and art.
- ◎ The cooler is made of seamless stainless steel tube with a high heat-exchange rate and big volume.
- ◎ This product adopts imported apparatus, which is highly sensitive and cut off the power when the tank is nearly out of water with a stable performance to protect the user.



### 【YNZD Series Technical Data 】

Specification	YNZD-3	YNZD-5	YNZD-10	YNZD-20
Output of water	3L/h	5L/h	10L/h	20L/h
Consumption Power	2.5kw	4.5kw	7.5kw	15kw
Mains voltage	220V	220V	380V	380V
Working size	33×24×62(cm)	39×28×76(cm)	42×30×87(cm)	50×40×83(cm)

### 【YNSD Series Technical Data (Automatic cut-off Water distiller)】

Specification	YNSD-3	YNSD-5	YNSD-10	YNSD-20
Output of water	3L/h	5L/h	10L/h	20L/h
Consumption Power	2.5kw	4.5kw	7.5kw	15kw
Mains voltage	220V	220V	380V	380V
Working size	33×24×62(cm)	39×28×76(cm)	42×30×87(cm)	50×40×83(cm)