



Performance Introduction

- Integrated aviation composite table top.
- Brand-new sheet metal housing with powder coating.
- The control system can control the lighting, the film viewer, the preheater of the laryngoscope and the patient chair.



坐墊顏色選擇 天藍 象牙黃 黑色

Standard Configuration

- | | |
|--|------------------------------|
| 1 Tabletop 1 PC | 8 Vacuum pump 1 PC |
| 2 Main body 1 SET | 9 Instrument tray 1 PC |
| 3 Spray gun 3 PCS (2 bent, 1 straight) | 10 Tweezer cup 1 PC |
| 4 Suction gun 1 PC | 11 Cotton cup 2 PCS |
| 5 Laryngoscope preheater 1 SET | 12 Sewage system 1 PC |
| 6 LED illuminating light 1 set | 13 Microcomputer system 1 PC |
| 7 Compressor 1 PC | 14 Doctor stool 1 PC |

Optional

- Filme Viewer 1PC
- Patient Chair 1PC



Colour: Blue Ivory

◆ Basic type patient chair technical Parameters

- 1.Voltage: AC220V+10% 50HZ
- 2.Power: 230W(max), low power consumption, energy saving.
- 3.Under weight pressure: 250kg(max)
- 4.Seat height adjustment: 480mm(minimum)- 610mm (maximum)
- 5.The backrest adjusting angle: 90°
- 6.Pillow flex: 100mm (foldable, detachable)

◆ Standard type patient chair technical Parameters

- 1.Voltage: AC220V 50HZ
- 2.Power: 230W(max), low power consumption, energy saving.
- 3.Under weight pressure: 250kg(max)
- 4.Seat height adjustment: 495mm(minimum)- 615mm (maximum)
- 5.The backrest adjusting angle: 90° -135°
- 6.Pillow flex: 100mm (foldable, detachable)
7. The seat can be rotated 180°



Colour: Blue Ivory



Colour: Blue Ivory

◆ Luxurious Type patient chair technical Parameters

- 1.Voltage: AC220V+10% 50HZ
- 2.Power: 280W(max), low power consumption, energy saving.
- 3.Under weight pressure: 250kg(max)
- 4.Seat height adjustment: 550 mm(minimum)- 630mm (maximum)
- 5.The backrest adjusting angle: 90° -170°
- 6.Pillow flex: 100mm (foldable, detachable)
7. The seat can be rotated 180°

