

Standard Configuration:

1,	Marble countertops	1 pcs	8. Lighting system	1 set
2,	Box	1 set	9. Compressor	1 Unit
3,	Spray gun		10. Vacuum pump	1 Unit
(2	straight and 1 bent)	1 handfuls	11. Instrument trays	2 pcs
4,	Sniffing guns	1 handful	12. Tweezer cups	2 pcs
5,	Blow guns	1 handful	13, Cotton ball cups	4 pcs
6.	Preheater	1 set	14. Drop bottles	2 pcs
7,	Key touch the key	1 set	15. Writing table	1 sheet

 16. Instrument recovery tray
 1 pcs

 17. Dirt buckets
 1 pcs

 18. Endoscope sterilization cup
 1 pcs

 19. Sewage discharge
 1 set

 20. Ionizer
 1 set

 21. Control system
 1 set

 22. Doctor's chair
 1 sheet

Options:

- High-definition endoscopic camera system
- PC (LCD monitor/dedicated endoscope software)
- Printer
- Light guides
- LED Cold light source
- Viewing lights

- Rigid tubes:
- ¢ 2.7/4mm,0°/30° Ear endoscopy ¢ 2.7/4mm,0°/30° Nasa1 endoscopy ¢ 8mm,70°/90° Laryngea1 endoscopy
- Electric patient chair
- LED headlight
- Built-in cold light source
- Disposable suction tubes (Φ2.0, Φ3.0)