



FETAL DOPPLER

Features

- Easy and convenient operation
- High-fidelity, crystal clear sound
- Earphone and speaker are feasible
- High sensitivity doppler probe
- Low ultrasound dosage
- Display with color LCD
- FHR with voice function(option)
- Bluetooth APP(option)

Fetal Doppler

Technical Specification

Ultrasound frequency : 2MHz/2.5MHz/3MHz

Ultrasound intensity : < 10W/cm²

Power supply:

2pcs 1.5V AA battery

2pcs DC Ni-Mh rechargeable battery

AC 220/110V.50/60Hz(for D and C)

FHR Measuring range: 50~230bpm

FHR resolution:1bpm

FHR accuracy: ±2bpm

Power consumption:< 1W

610 Series



Function	610B	610D	610G	610C
Display	B/W	B/W	Color	Color
Backlight	Yes	Yes	Yes	Yes
Build-in speaker	Yes	Yes	Yes	Yes
Various work modes	No	No	Yes	Yes
Volume adjustment	Yes	Yes	Yes	Yes
Battery indicator	Yes	Yes	Yes	Yes
FHR display	Yes	Yes	Yes	Yes
FHR curve	No	No	Yes	Yes
FHR voice function	Option	Option	Option	Option
Alkalinity battery	Yes	No	Yes	No
Ni-Mh battery	No	Yes	No	Yes
Charge in device	No	Yes	No	Yes
Earphone	Option	Option	Option	Option
Carry bag	Option	Option	Option	Option
Bluetooth	Option	Option	Option	Option

All specifications are subject to change without notice.