

## Oil Bath

### HH-4

- ◆ Microprocessor PID control
- ◆ Digital display setting and actual temperature
- ◆ Digital display running timer
- ◆ Annular butterfly cover, easy to operate
- ◆ High efficiency heat preservation material, small heat loss

### Safety Mechanism:

- ◆ PID self-tuning , over temperature alarm
- ◆ Automatic alarm after operation
- ◆ Electronic timing, PID self - setting, overtemperature alarm
- ◆ Automatic alarm after work
- ◆ Temperature control protection and alarm function



HH-4

Model	HH-1	HH-2	HH-4	HH-6	HH-8
Capacity	3L	6L	10L	16L	22L
Working size(mm)	160x160x120	300x160x120	300x300x120	460x300x120	600x300x120
Power	500W	1000W	1500W	2000W	2500W
Temperature fluctuation	±0.5°C	±0.5°C	±1°C	±1.5°C	±1.5°C
Temperature range			RT+5-300°C		
Temperature homogeneity			±1°C		
Voltage	110V 60Hz /220V 50Hz				
Dimension(mm)	330x230x240	470x230x240	470x370x240	630x370x240	770x370x240
Package(mm)	410x280x290	530x290x280	530x430x280	690x430x280	840x450x280
G.W./N.W.	6/5.5kg	9/8kg	11/10kg	14/13kg	17/16kg

## Oil Bath

### HH-Y8

- ◆ High precision microcomputer to control temperature , accurate and reliable temperature control
- ◆ The set temperature and actual temperature are displayed in two positions
- ◆ Internal and external circulation, internal circulation to ensure uniform and constant temperature
- ◆ High efficiency insulation material, low heat loss
- ◆ With higher and lower temperature alarm setting functio



HH-Y8

Product name	Soil oil bath	Oil bath	Oil bath	Cycle oil bath
Model	HH-Y8	HH-240	HH-280	602
Working size (mm)	600x300x120	240x240x240	280x280x300	400x280x210
Power	2500W	1500W	2000W	2000W
Temperature range	RT+5-300°C	RT+5-300°C	RT+5-300°C	RT+5-300°C
Temperature fluctuation	±1.5°C	±1°C	±1°C	±0.5°C
Temperature accuracy	±1.5°C	±1°C	±1°C	±0.5°C
Voltage	110V 60Hz /220V 50Hz			
Dimension (mm)	770x370x360	320x320x360	360x360x420	480x450x430
Package (mm)	840x430x480	380x380x420	420x420x480	540x510x490
G.W./N.W.	21/19kg	12/10kg	14/12kg	20/18kg