













IB-9052A

- ♦ Microprocessor control
- ◆Digital display setting and actual temperature
- ♦ Over temperature alarm, PID self-tuning
- ♦ Digital display running timer,
- ◆Stainless steel interior
- ♦ With glass inner door, With hidden door locks

Safety Mechanism:

- ◆Electricity leakage protection
- ◆ Mechanical overheat protection, cut off the heating sources
- ◆ Automatic alarm after operation

Scope of application:

This product is used for bacterial culture, breeding, fermentation and other constant temperature experiments in scientific research departments such as medical and health, pharmaceutical industry, biochemistry, industrial production and agricultural science.



Model	IB-9025A	IB-9052A	IB-9082A	IB-9162A	IB-9272A					
Capacity	15L	50L	80L	160L	270L					
Working size (mm)	250x250x250	350x350x410	400×400×500	500x500x650	600x600x750					
Power	150W	250W	300W	400W	600W					
Temperature range	Room temperature+5~65℃									
Temperature resolution	0.1℃									
Temperature fluctuation	±0.5℃									
Timer	0-999min									
Voltage	110V 60Hz /220V 50Hz									
Dimension(mm)	380x370x580	480x470x740	530x530x840	640x630x1000	730x720x1090					
Package(mm)	490x470x640	580x580x800	630x650x980	730×710×1120	880x860x1280					
G.W./N.W.	35/21kg	46/35kg	62/47kg	90/64kg	112/86kg					

PIOWAY

Laboratory Incubator

DNP-9052A

- ◆Microprocessor control
- ◆Digital display setting and actual temperature
- ◆Over temperature alarm, PID selftuning
- ◆Digital display running timer,
- ◆Stainless steel interior
- ◆With glass observe window

Safety Mechanism:

◆Automatic alarm after operation

Scope of application:

This product is used for bacterial culture, breeding, fermentation and other constant temperature experiments in scientific research departments such as medical and health, pharmaceutical industry, biochemistry, industrial production and agricultural science.



Model	DNP-9025A	DNP-9022A	DNP-9030A	DNP-9052A	DNP-9082A	DNP-9162A	DNP-9272A				
Capacity	15L	25L	30L	50L	80L	160L	270L				
Working	250x250	300x300	300x300	350x350	400x400	500x500	600x600				
size (mm)	x250	x280	x350	x410	x500	x650	x750				
Power	150W	180W	200W	250W	300W	400W	600W				
Temperature range	Room temperature+5~65°C										
Temperature resolution				0.1℃							
Temperature fluctuation	±0.5℃										
Timer	0-999min										
Voltage	110V 60Hz /220V 50Hz										
Dimension	370x380	420x430	420x430	470x490	520x540	620x640	720x740				
(mm)	x540	x570	x640	x700	x800	×940	x1010				
Package	460x470	510x520	510x520	570x580	620x630	720x730	810x820				
(mm)	x610	x640	x720	x770	x900	x1080	×1170				
G.W./N.W.	30/20kg	35/25kg	38/28kg	48/36kg	59/44kg	77/60kg	103/82kg				