



# Sterilizer



## Hot-air Sterilizer

### GRX-9053A

- ◆ Microprocessor control
- ◆ With double glass window
- ◆ Digital display setting and actual temperature
- ◆ Digital display running timer
- ◆ PID self-tuning, over temperature alarm
- ◆ Automatic alarm after operation
- ◆ Stainless-steel interior

### Ideal For:

- ◆ Hot air sterilizer uses high temperature dry heat to oxidize microorganisms, protein denaturation, and poisoning caused by dielectric concentration.
- ◆ Destroy cell protoplasm by oxidation, make microorganisms die, and kill all microorganisms within a certain heating time



Model	GRX-9015A	GRX-9023A	GRX-9053A	GRX-9073A	GRX-9123A	GRX-9203A
Capacity	15L	25L	50L	70L	105L	225L
Working size(mm)	250x250x250	300x300x275	415x375x345	450x350x450	550x370x520	600x500x750
Power	300W	500W	1000W	1200W	1500W	2500W
Temperature range	RT+10-250°C					
Temperature resolution	0.1°C					
Temperature fluctuation	±1°C					
Timer	0-999min					
Voltage	110V 60Hz / 220V 50Hz					
Dimension(mm)	520x450x440	570x530x460	690x610x530	610x700x620	720x700x720	840x820x830
Package(mm)	600x480x520	660x600x530	790x660x620	720x780x730	820x750x810	940x870x1130
G.W./N.W.	33/22kg	41/28kg	50/39kg	62/48kg	78/58kg	114/85kg