

## Non-Vacuum Blood Collection Tube

- ★ Made of medical grade PP material.★ Siliconization treatment on the inner surface, no blood adhesion after collection.
- ★ High quality gel makes the serum separation complete with no fibrin adhering to the inner surface of the tube.
- There uniformLy sprays the quantitative coagulant/anticoagulant, which will mix smoothly and completely with the blood, thus ensures fast coagulation or effective anticoagualtion. 

  Good compatibility makes the tube match with various kinds of biochemical, immunity and haematology analyzers.
- \* Labes which can indicate the precise blood collection volume can be customized according to clients requirement.











REF	Original REF.	Tube Specification	Nominal Volume	Additive	Qty/Cs
0302-0212	— КЈ203 —	Ø13×75mm	2mL	EDTA.K2	100PCS×20
0302-0522		Ø13×100mm	5mL	EDTA.K2	100PCS×20
0307-0312	KJ204	Ø13×75mm	3mL	Lithium Heparin	100PCS×20
0307-0522		Ø13×100mm	5mL	Lithium Heparin	100PCS×20
0306-0312		Ø13×75mm	3mL	Sodium Heparin	100PCS×20
0306-0522		Ø13×100mm	5mL	Sodium Heparin	100PCS×20
0301-0312	— KJ204-1 —	Ø13×75mm	3mL	Gel & Clot Activator	100PCS×20
0301-0522		Ø13×100mm	5mL	Gel & Clot Activator	100PCS×20
0308-0212	— KJ204-2 —	Ø13×75mm	2mL	Sodium Fluoride+Heparin Sodium	100PCS×20
0308-0522		Ø13×100mm	5mL	Sodium Fluoride+Heparin Sodium	100PCS×20
0314-0212		Ø13×75mm	2mL	Sodium Fluoride+EDTA.K2	100PCS×20
0314-0522		Ø13×100mm	5mL	Sodium Fluoride+EDTA.K2	100PCS×20
0304-0312	— КЈ204-4 —	Ø13×75mm	3mL	No Additive	100PCS×20
0304-0522		Ø13×100mm	5mL	No Additive	100PCS×20

- \* KANGJIAN now lauch the new type non-vacuum tube.
- \* With the excellent design for the stopper and tube structure, it is convenient for inserting or removing the stopper.
- ★ Perfect match between the stopper and the tube ensures no leakage.
- \* Easy and convenient for one-hand operation during the experiment.



