



Plasma Thawing Instrument

RJ-8

- ◆ Microprocessor PID control
- ◆ Digital display actual temperature
- ◆ Stainless steel interior
- ◆ Stainless steel basket
- ◆ Automatic arms swing system
- ◆ Optional stainless steel flat cover or roof cover

Safety Mechanism:

- ◆ With over temperature protection, alarm and auto cut-off
- ◆ Prevent dry burning heating components, effectively protect heating elements



Model	RJ-2	RJ-4	RJ-8	RJ-12
Capacity	2 bags	4 bags	8 bags	12 bags
Heating power	800W	1000W	1500W	1800W
Temperature range		37°C ± 1°C		
Temperature fluctuation		± 0.5°C		
Defrosting time		≤ 10min		
Shaking frequency		80rpm		
Voltage	110V 60Hz / 220V 50Hz			
Dimension(mm)	450x340x410	450x460x410	690x460x410	690x550x410
Package(mm)	600x420x490	600x560x490	690x560x490	840x650x490
G.W./N.W.	22/16kg	29/21kg	41/30kg	61/46kg