













## Water Bath Shaker

## SHA-C

- ◆Microprocessor PID control
- ◆Digital display actual temperature
- ♦With mechanical timer
- ◆Good insulation performance with stainless steel roof cover
- ◆Stainless steel interior and cover
- ♦Universal spring clamp, use for all kinds of reagent bottle
- ◆Can optional speed and time digital display control(SHZ-88A)

## Scope of application:

◆ Mainly applicable to scientific research departments such as universities and colleges, medical, petrochemical, health and epidemic prevention, environmental monitoring, etc.



Model	SHA-C	SHZ-82	SHA-B	SHZ-88A
Shaking mode	Reciprocating	Orbital	Reciprocating & Orbital	Reciprocating
Working size(mm)	490x390x160			490x390x160
Clamp size(mm)	390x290			390x290
Temperature	RT+5~99℃			RT+5~99°C
Temperature fluctuation	±0.5℃			±0.2℃
Heating power	2000W			2000W
Timer range	0-999min or with no limit			0-999min
Shaking frequency	Start up -300rpm, stepless speed adjustment			40-280rpm
Amplitude	20mm			26mm
Voltage	110V 60Hz /220V 50Hz			
Dimension(mm)	710x510x450	710x510x450	710x510x540	680x500x400
Package(mm)	800x600x550	800x600x550	800x600x620	770×590×480
G.W./N.W.	50/35kg	50/35kg	55/40kg	50/35kg