

## Cycle Water Bath

### ZEN-TW

- ◆ Microprocessor control, touch panel control
- ◆ Multi purpose: room temperature to 99°C
- ◆ Optional function: room temperature to 200 °C
- ◆ One-timer-molding cistern, easy to clean
- ◆ LCD display temperature setpoint and temperature measurements, with 0.1 °C resolution
- ◆ Inside and outside pump cycle, inside circulation ensures temperature uniformity



ZEN-TW

Model	ZEN-TW
Capacity	27L
Inner size(mm)	460x300x200
Working size(mm)	340x300x200
Power	1200W
Temperature range	RT+5~99°C
Temperature accuracy	±0.5°C
Pump circulation mode	Inner and outer cycle
Pump flow	8-10L/min
Voltage	110V 60Hz / 220V 50Hz
Dimension(mm)	530x330x430
Package(mm)	600x390x490
G.W./N.W.	16/14kg

### Safety Mechanism:

- ◆ Automatic alarm after work
- ◆ With higher and lower temperature alarm setting function
- ◆ Temperature control protection and alarm function

## Cycle Water Bath

### WH-4C/WH-6C/WH-60C

- ◆ Microprocessor control
- ◆ Digital display inner temperature
- ◆ One-timer-molding cistern
- ◆ Good insulation performance with stainless steel cover
- ◆ With inner cycle
- ◆ Optional stainless steel flat cover or roof cover



Model	WH-4C	WH-6C	WH-60C
Capacity	12L	16L	62L
Power	1000W	1500W	2500W
Pump power	40W	40W	60W
Circulation pump flow	≥4L/min	≥4L/min	≥8-10L/min
Temperature fluctuation	±0.1°C	±0.1°C	±0.5°C
Temperature range	RT+5 -99°C	RT+5 -99°C	RT+5 -99°C
Voltage	110V 60Hz /220V 50Hz	110V 60Hz /220V 50Hz	110V 60Hz /220V 50Hz
Working size(mm)	320x320x120	460x300x120	620x500x200
Dimension(mm)	460x330x240	640x330x250	700x450x320
Package(mm)	530x390x280	710x390x290	800x550x370
G.W./N.W.	9/8kg	12/11kg	49/35kg