



Water Bath

HH-S2/HH-S4/HH-S6

- ◆ Microprocessor PID control
- ◆ Digital display setting and actual temperature
- ◆ Annular butterfly cover, easy to operate
- ◆ High efficiency heat preservation material, small heat loss



Model	HH-S1	HH-S2	HH-S4	HH-S6	HH-S8
Capacity	3L	6L	10L	16L	22L
Power	300W	500W	1000W	1500W	2000W
Temperature fluctuation	±0.5°C	±0.5°C	±1°C	±1.5°C	±1.5°C
Temperature range			RT+5 -99°C		
Temperature homogeneity			±1°C		
Voltage			110V 60Hz /220V 50Hz		
Working size(mm)	160x160x120	300x160x120	300x300x120	460x300x120	600x300x120
Dimension(mm)	230x230x240	370x230x240	370x370x240	530x370x240	670x370x240
Package(mm)	290x290x280	410x280x280	420x410x280	610x430x280	730x430x280
G.W./N.W.	3.5/3kg	4.5/4kg	6/5 kg	7.5/6.5kg	9/8kg

Water Bath

WH-1/WH-2/WH-4/WH-6

- ◆ Microprocessor control,digital display inner temperature
- ◆ One-timer-molding cistern
- ◆ Stainless steel roof cover (optional)

Scope of application:

- ◆ Suitable for thermostatic , baking, dewaxing, sterilization



Model	WH-1	WH-2	WH-4	WH-6
Capacity	3L	6L	12L	16L
Power	200W	500W	1000W	1500W
Temperature fluctuation	±0.5°C	±0.5°C	±1°C	±1.5°C
Temperature range	RT+5 -99°C			
Temperature homogeneity		±1°C		
Voltage		110V 60Hz /220V 50Hz		
Working size(mm)	150x140x120	300x150x120	320x320x120	460x300x120
Dimension(mm)	180x180x240	330x180x240	360x330x240	540x330x250
Package(mm)	230x220x300	380x230x280	410x400x280	590x400x310
G.W./N.W.	3/2.5kg	4/3 kg	6/5 kg	8/7 kg