

DIRECT OPHTHALMOSCOPE

YZ11  

Direct Ophthalmoscope

- Elegant appearance.
- It is convenient for mobile inspection.
- Both dry battery and rechargeable battery power are available.



Specifications

Illumination Type	Large Spot , Middle Spot , Small Spot , Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±15 ±20 -25 -35D
Illumination Source	2.6V/0.78W , Incandescent Bulb
Power Consumption	C Alkaescent battery , rechargeable battery

YZ11D  

Rechargeable Direct Ophthalmoscope

- Improved optical design and excellent optical performance.
- Imported halogen bulb provides high and adjustable brightness and good color rendering.
- Smart, portable and convenient for mobile inspection.
- Equipped with special recharger quickly and conveniently recharging, against over recharging.
- Anti-recharging and anti-discharging function protects the battery.



Specifications

Illumination Type	Large Spot , Small Spot , Slit , Net , Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±15 ±20 -25 -35D
Illumination Source	3.5V/2.8W Halogen Bulb
Power Supply	AC110V/220V±10% , 50Hz/60Hz
Power (Main Body of Ophthalmoscope)	Li Batteny Pack 3.7V , Rechargeable